# INFRARED

## BODY THERMOMETER

**NON-CONTACT**

PITHG01

USER MANUAL

## EN User manual

### Non-contact electronic infrared body thermometer

This manual is part of the product and should be stored in a way that allows you to read it in the future. It contains important information about the settings and operation of the device.

The set includes:

- Electronic thermometer equipped with an infrared sensor
- User manual

Precautions:

Before using the thermometer for the first time, please read the following notes:

- This product is only intended for the use described in this manual
- Measurement performed with this product cannot substitute for a proper medical examination
- The infrared sensor is the most delicate part of the thermometer. Do not touch the lens of the infrared sensor with your fingers.
- If the product is not used for more than two months, remove the batteries and store them separately.
- Dispose of used batteries in accordance with applicable law. Batteries may explode if overheated. After discharged, they must be collected separately from other waste.
- Do not submerge the product in water or expose it to sunlight, which may damage it.
- Do not drop or expose to electrical surges.
- Do not use outdoors or near a fireplace or stove.
- If necessary, clean the surface of the glass with the infrared sensor with 70% isopropyl alcohol.
- The device should be used in an ambient temperature between 5 ° C and 40 ° C (0 ° F to 140 ° F). The thermometer should not be used outside this range. To get the most accurate measurement results, the thermometer and the test person should be in the same room for about 15-20 minutes to equilibrate their temperature before the measurement.
- Maintenance tasks that you can perform only include battery replacement and cleaning.

Thermometer application:

This product measures body temperature by measuring the infrared radiation emitted from the skin surface on the forehead. The non-contact thermometer is designed to measure body temperature (Body mode) and surface temperature (Surface mode). It can also check the temperature of the room, food, bath water, baby milk, etc.

Product description and method of operation:

The thermometer consists of a plastic housing, LCD display, function buttons, infrared temperature sensor, components and electronic circuits.

The non-contact thermometer is a handheld, battery-powered device designed to periodically measure and monitor the temperature of the human body by measuring the infrared radiation emitted from the skin surface on the forehead of infants, children and adults. The optical system of the thermometer directs and focuses the infrared energy to the detector, and then the electronics of the instrument convert the signal from the detector to the reading shown on the display.


Drawing and description of the device structure:

1. Infrared sensor
2. LCD display
3. DOWN button
4. UP button
5. SET button
6. ON / OFF button - trigger
7. Battery socket

Display drawing and description:

1. Surface mode - measuring surface temperature
2. Body mode - measuring body temperature
3. Measurement result
4. Measurement result between 32 ° C - 37.2 ° C - ☺ icon
5. Measurement result over 37.2 ° C - ☹ icon
6. Low battery indicator
7. Icon of speaker and volume level
8. Setting the scale of the measurement result to ° C or ° F
9. Number of the saved measurement

Taking measurements


Comments:

- Measurements should be made only on the exposed part of the forehead, any covering (i.e. hair or a hat) may cause an erroneous measurement.
- Accurate measurement requires cleaning the test surface (i.e. forehead) from make-up and other substances that could distort the measurement result. Bring any hair off your forehead, wipe away any sweat.

- If you are testing the temperature of babies who are usually wrapped in a blanket or clothes, the neck should be exposed for at least 1 minute before measuring.
- Cold compresses and other methods of reducing fever may affect the correct course of the measurement and its results.
- Measurements should be made in an environment with a stable temperature distribution. Do not measure near an operating fan or air conditioner, in a strong blast of cold air.
- Do not take measurements in the area of strong sunlight.
- Hold the thermometer approximately 3-5 cm from your forehead, press the ON / OFF scan button. The temperature will be displayed immediately.
- If the measurement result seems to you inadequate to the actual body temperature you feel, repeat the measurement with the sensor pointing at the skin behind the ear.

The measuring procedure:

1. Aim the thermometer at the forehead, squeeze the trigger (ON / OFF button). The LCD will temporarily show all the figures to calibrate the ambient temperature during this time.
2. After completing the initial calibration process, the device will emit a sound signal and the LCD screen will display the measurement result, simultaneously assuming one of the three colors (green, yellow, red) according to the saved settings regarding the range of the measurement result.

Changing the measurement mode:

1. Turn on the device. Press the "SET" button to toggle between Body mode and Surface temp. Set to "Surface" mode.
2. The thermometer should be brought closer to the surface (i.e. liquid) at a distance of 5-15 cm. Do not immerse the thermometer.
3. The measuring range is 0 - 100 ° C.

## Thermometer calibration

It is possible to calibrate the thermometer on the basis of indications of another reference thermometer.

1. Calculate the difference of readings between the reference thermometer and the Platinet thermometer.
2. Turn on the device, press the "SET" button and hold for 3 seconds, the LCD will show F1.
3. 0.0 flashes on the display. Use the + or - buttons to set the difference value. Each time the buttons are pressed, the value will increase / decrease by 0.1 ° C.
4. Confirm the settings by pressing the SET button

## High temperature alert display

The high temperature alert display apply only to the Body mode, i.e. body temperature measurement. When taking surface measurements in Surface mode, each measurement will have a green LCD backlight. To set the alarm temperature:

1. Turn on the instrument, press the "SET" button and hold for 3 seconds, the LCD shows F1, press the "SET" button again and go to F2.
2. Use the "UP" and "DOWN" buttons to set the temperature level, when exceeded, an alarm will be triggered and the LED screen will illuminate red.

## Switching the scale of the displayed measurements from ° C to ° F

1. Turn on the instrument, press the "SET" button and hold for 3 seconds, the LCD shows F1, press the "SET" button two more times and go to F3.
2. Use the "UP" and "DOWN" buttons to switch the scale unit in which you want to read the temperature measurement results.

## Enabling and disabling the acoustic signal





1. Turn on the instrument, press the "SET" button and hold for 3 seconds, the LCD shows F1, press the "SET" button three more times and go to F4.
2. Use the "UP" and "DOWN" buttons to turn the sound signal on or off.

## Thermometer memory:

Temperature measurement results are saved automatically, the thermometer can store up to 32 values, and they will be overwritten with new data when the thermometer memory is full.

To browse the device memory, just turn it on and press the "UP" or "DOWN" buttons. If there are no values saved in the device memory, the symbol "---" will appear.

The meaning of the symbols that may appear on the display:

	The battery is running low, replace the batteries with new ones
	EEPROM error, send the device back for servicing
	Infrared sensor error, send the device for servicing
	Temperature measurement above 42.2 ° C
	Temperature measurement result below 32 ° C
	Ambient temperature is too high for measurement
	Ambient temperature is too low for measurement

Cleaning and storage:

- After use, wipe the thermal sensor with a damp cotton swab, dry with a paper towel and store in the package.
- The device should be cleaned once a week with a cloth soaked in water or a neutral detergent, avoiding to flood it or excessively wet it.
- The device should not be repaired by yourself. Do not modify this thermometer without the approval of the manufacturer. Do not disassemble or attempt to repair the device or its component yourself.

Battery change

The device requires two AAA (LR03) batteries. Open the battery compartment. To do this, slide down the battery compartment cover on the back. Take out the used batteries from the battery compartment. Install new batteries. Make sure the batteries are inserted the right way round. Close the battery compartment.

Security

Temperature measurement results should not be used as the basis for medical examination. For measuring babies, keep the baby still. Movement can cause erroneous readings. It is not allowed to modify or disassemble the device. If the thermometer is not used for an extended period of time, remove the batteries and store in a safe place. The device is not suitable for use in the presence of flammable anesthetics or oxygen. This device has been tested and found to comply with the limits for a medical device in EN 60601-1: 2006 + A11: 2011 + A1: 2013 + A12: 2014. The limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to other equipment in the vicinity.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause interference to other equipment, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Change the direction or position of the receiving device.
- Increase the separation between the devices.
- Connect the equipment into an outlet on a circuit different from that to which the other equipment or devices are connected.
- Consult the manufacturer's service technician for assistance. Do not measure during or immediately after breastfeeding. Temperature readings may vary by location. Therefore, we recommend that you measure the temperature at 2 points. Do not use the thermometer in environments with high humidity.

Technical parameters:

Product name: non-contact infrared thermometer

Operating ambient temperature: 16 °C -35 °C

Working atmospheric pressure: 70kPa-106kPa

Working air humidity: ≤85%

Supply voltage: 3.0 V DC (2 AAA alkaline batteries)

Size: 149 \* 77 \* 43mm

Weight: 105g / 131g

Measuring range: Body temperature: 32 °C - 42.9 °C

Measurement accuracy: 35.0-42.0 °C: ± 0.2 °C; outside the range of 35.0-42.0 °C: ± 0.3 °C

Measurement distance: 3-5 cm

Measurement time: <2 seconds

3-color display: body temperature: 34.0-37.3 °C green; 37.4-38 °C yellow; 38.1-42.6 °C red

Automatic shutdown: 15 seconds

Storage and transportation: -20 °C -55 °C, Relative Humidity: ≤93%

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